Using *crlmm* to genotype data from Illumina's Infinium BeadChips

Matt Ritchie

October 7, 2010

1 Getting started

In this user guide we read in and genotype data from 40 HapMap samples which have been analyzed using Illumina's 370k Duo BeadChips. This data is available in the hapmap370k package. Additional chip-specific model parameters and basic SNP annotation information used by CRLMM is stored in the human370v1cCrlmm package. The required packages can be installed in the usual way using the biocLite function.

```
> source("http://www.bioconductor.org/biocLite.R")
> biocLite(c("crlmm", "hapmap370k", "human370v1cCrlmm"))
```

2 Reading in data

The function readIdatFiles extracts the Red and Green intensities from the binary idat files output by Illumina's scanning device. The file samples370k.csv contains information about each sample.

Reading in this data takes approximately 90 seconds and peak memory usage was 1.2 GB of RAM on our linux system. If memory is limiting, load the ff package and run the same command. When this package is available, the objects are stored using disk rather then RAM. The RG object is an NChannelSet which stores the Red and Green intensities for each bead-type. The scanning date of each array is stored in protocolData.

```
> class(RG)
[1] "NChannelSet"
attr(,"package")
[1] "Biobase"
> dim(RG)
Features
          Samples
  381079
                40
> slotNames(RG)
[1] "assayData"
                          "phenoData"
[3] "featureData"
                          "experimentData"
[5] "annotation"
                          "protocolData"
[7] ".__classVersion__"
> channelNames(RG)
[1] "G"
            "R"
                   "zero"
> exprs(channel(RG, "R"))[1:5, 1:5]
      4030186347_A 4030186263_B 4019585415_B
10008
                321
                              170
                                           2961
               1738
10010
                             3702
                                           3105
10025
                 80
                              101
                                            145
10026
               5043
                             1856
                                           6519
10039
               4905
                             2464
                                           9080
      4031058127_B 4031058211_B
10008
               3468
                              262
10010
               3425
                               70
10025
                 29
                               21
10026
               8304
                             9872
10039
               9788
                            10867
```

> exprs(channel(RG, "G"))[1:5, 1:5]


```
10008
                            4484
              4183
                                         3765
10010
              2593
                              51
                                         3824
              2768
                            2322
                                         3435
10025
10026
               216
                            2840
                                          211
10039
               297
                                          345
                            3016
      4031058127_B 4031058211_B
10008
              3558
                            6502
              3528
10010
                            6154
10025
                            3608
              3471
10026
               164
                             188
10039
               361
                             380
> pd = pData(RG)
> pd[1:5, ]
             HapMap.Name Gender
                                         Plate
                 NAO6991 Female WG1000442-DNA
4030186347_A
4030186263_B
                 NAO7000 Female WG1000442-DNA
                 NA10859 Female WG1000453-DNA
4019585415_B
4031058127_B
                 NA11882 Female WG1000453-DNA
4031058211_B
                 NA06993
                           Male WG1000447-DNA
             Well SentrixPosition
4030186347_A E11
                     4030186347_A
4030186263_B D08
                     4030186263_B
4019585415_B
              B02
                     4019585415_B
4031058127_B D08
                     4031058127_B
4031058211_B D11
                     4031058211_B
> scandatetime = strptime(protocolData(RG)[["ScanDate"]],
      "%m/%d/%Y %H:%M:%S %p")
> datescanned = substr(scandatetime, 1,
      10)
> scanbatch = factor(datescanned)
> levels(scanbatch) = 1:16
> scanbatch = as.numeric(scanbatch)
```

4030186347_A 4030186263_B 4019585415_B

Plots of the summarised data can be easily generated to check for arrays with poor signal.

```
> par(mfrow = c(2, 1), mai = c(0.4, 0.4, 0.4, 0.4), oma = c(1, 1, 0, 0))
> boxplot(log2(exprs(channel(RG, "R"))),
```

```
+ xlab = "Array", ylab = "", names = 1:40,
+ main = "Red channel", outline = FALSE,
+ las = 2)
> boxplot(log2(exprs(channel(RG, "G"))),
+ xlab = "Array", ylab = "", names = 1:40,
+ main = "Green channel", outline = FALSE,
+ las = 2)
> mtext(expression(log[2](intensity)), side = 2,
+ outer = TRUE)
> mtext("Array", side = 1, outer = TRUE)
```



3 Genotyping

Next we use the function **crlmmIllumina** which performs preprocessing followed by genotyping using the CRLMM algorithm.

> dim(crlmmResult)

[1] "SnpSet"
attr(,"package")
[1] "Biobase"

Features Samples 346451 40

> slotNames(crlmmResult)


> calls(crlmmResult)[1:10, 1:5]

```
1 2 3 4 5 rs12354060 1 1 3 3 3 3 rs6650104 1 1 1 1 1 1 1 rs12184279 1 1 1 1 1 1 rs3115860 2 1 1 2 2 rs3115850 1 2 2 1 1 rs7515489 3 3 1 1 1 1 rs12124819 1 2 2 1 1 rs17160939 1 1 1 1 1 rs12086311 3 3 3 3 3
```

Plotting the SNR reveals no obvious batch effects in this data set (different symbols are used for arrays scanned on different days).

```
> plot(crlmmResult[["SNR"]], pch = scanbatch,
+ xlab = "Array", ylab = "SNR", main = "Signal-to-noise ratio per array",
+ las = 2)
```

Signal-to-noise ratio per array



4 System information

This analysis was carried out on a linux machine with 32GB of RAM using the following packages:

> sessionInfo()

R version 2.12.0 alpha (2010-09-21 r52960) Platform: x86_64-unknown-linux-gnu (64-bit)

locale:

- [1] LC_CTYPE=en_US.iso885915
- [2] LC_NUMERIC=C
- [3] LC_TIME=en_US.iso885915
- [4] LC_COLLATE=en_US.iso885915
- [5] LC_MONETARY=C
- [6] LC_MESSAGES=en_US.iso885915
- [7] LC_PAPER=en_US.iso885915
- [8] LC_NAME=C
- [9] LC_ADDRESS=C
- [10] LC_TELEPHONE=C
- [11] LC_MEASUREMENT=en_US.iso885915

[12] LC_IDENTIFICATION=C

attached base packages:

- [1] tools stats graphics grDevices
- [5] utils datasets methods base

other attached packages:

- [1] human370v1cCrlmm_1.0.1 hapmap370k_1.0.0
- [3] crlmm_1.7.15 oligoClasses_1.11.8
- [5] Biobase_2.9.1 weaver_1.15.0
- [7] codetools_0.2-2 digest_0.4.2

loaded via a namespace (and not attached):

- [1] affyio_1.17.4 annotate_1.27.1
- [3] AnnotationDbi_1.11.5 Biostrings_2.17.47
- [5] bit_1.1-4 DBI_0.2-5
- [7] ellipse_0.3-5 ff_2.1-2
- [9] genefilter_1.31.2 IRanges_1.7.34
- [11] mvtnorm_0.9-92 preprocessCore_1.11.0
- [13] RSQLite_0.9-2 splines_2.12.0
- [15] survival_2.35-8 xtable_1.5-6